

BARIATRIC SLYDE™



belluscura™



- Folds flat and stacks for easy storage
- Restraints with buckles to secure the occupant for transport
- Corrugated polypropylene, double-wall construction is lightweight and durable
- Strong web handles at corners and sides provide multi-point control for two to seven carriers
- The length of the web handle allows an ergonomically correct upright posture for operators during stair descent

800LB
Capacity

1 YEAR
Warranty

The Belluscura Bariatric Evacuation Slyde provides a safe, efficient way to evacuate non-ambulatory patients from any multistory building using the stairways. It folds up and stores in stacks up to five units in a storage sleeve. The sleeve can be wall-mounted for quick access in key locations.



BARIATRIC SLYDE™



belluscura™

Features and Benefits

- Folds flat and stacks 5 deep for easy storage
- Easy set-up and instruction pictograms on each unit
- Four strong web handles for multipoint control by 2 to 7 operators
- Size permits continued stairwell traffic flow during occupant evacuation
- Functions as post-evacuation surge capacity device
- Lightweight for easier transports
- Buckle restraints are quick and easy to use¹
- Compatible with Belluscura Stairwell Belay System for heavy occupant evacuations.
- Wall-mounted storage options are space-saving and convenient

Accessories



Belay Box and Bolt
6400-030-000



Provides stable wall anchor



Stairwell Belay
6400-020-000



5 COUNT
6400-010-001

Slyde Sleeve

Specifications

Model Number	6410-000-000
Material	Corrugated polypropylene
Height	8 in. (20.3 cm)
Width²	24 in. (61 cm)
Length	85 in. (215.9 cm)
Weight	9 lb. (4.1 kg)
Folded Height	2 in. (5.1 cm)
Max. Weight Capacity^{3,4}	800 lb. (363 kg)

1. Passenger restraint set includes waist and chest restraints.
2. Dimensions are measured from the outermost edges of the main frame. Specifications are rounded to the nearest whole number. Conversions are calculated before rounding.
3. Operators must consider the weight of the passenger, equipment and accessories when determining the total weight on the product.
4. Maximum weight capacity is total weight distributed in accordance with basic human anatomy. Belluscura reserves the right to change specifications without notice.

ORDER INFORMATION

email: customerservice@belluscura.com
tel: +1 786-703-7996
toll free: +1 844-883-0940

Belluscura, LLC.
15455 Dallas Pkwy., Suite 600
Addison, TX 75001
United States
tel: +1 786-703-7996
toll free: +1 844-883-0940

Belluscura, Ltd.
Oxford Center for Innovation
New Road, Oxford OX1 1BY
United Kingdom
tel: +44 (0) 1865 261506
toll free: +1 844-883-0940