



The Belluscura Evacuation Slyde provides a safe, efficient way to evacuate non-ambulatory occupants from any multistory building using the stairways. It folds flat and stores in stacks of up to five units in a storage sleeve. The sleeve can be wall-mounted for quick access in key locations.



## belluscura<sup>™</sup>

- Folds flat and stacks for easy storage
- Restraints with buckles to secure the occupant for transport
- Corrugated polypropylene, double-wall construction is lightweight and durable
- Strong web handles at corners and sides provide multi-point control for two to seven carriers
- The length of the web handle allows an ergonomically correct upright posture for operators during stair descent

**500**LB Capacity







### **Features and Benefits**

- Folds flat and stacks 5 deep for easy storage
- Easy set-up and instruction pictograms on each unit
- Four strong web handles for multipoint control by 2 to 7 operators
- Size permits continued stairwell traffic flow during occupant evacuation

- Functions as post-evacuation surge capacity device
- · Lightweight for easier transports
- Buckle restraints are quick and easy to use<sup>1</sup>
- Compatible with Belluscura Stairwell Belay System for heavy occupant evacuations.
- Wall-mounted storage options are space-saving and convenient

### **Accessories**



### **Specifications**

**Model Number** 6400-000-000 **Material** Corrugated polypropylene 8 in. (20.3 cm) Height Width2 19 in. (48.3 cm) Length 85 in. (215.9 cm) Weight 9 lb. (4.1 kg) **Folded Height** 2 in. (5.1 cm) Max. Weight Capacity3,4 500 lb. (227 kg)



# belluscura"

- 1. Passenger restraint set includes waist and chest restraints.
- Dimensions are measured from the outermost edges of the main frame.
  Specifications are rounded to the nearest whole number. Conversions are calculated before rounding.
- 3. Operators must consider the weight of the passenger, equipment and accessories when determining the total weight on the product.
- 4. Maximum weight capacity is total weight distributed in accordance with basic human anatomy. Belluscura reserves the right to change specifications without notice.

## **ORDER INFORMATION**

email: customerservice@belluscura.com tel: +1 786-703-7996

toll free: +1 844-883-0940

#### Belluscura, LLC.

15455 Dallas Pkwy., Suite 600 Addison, TX 75001 United States

tel: +1 786-703-7996

toll free: +1 844-883-0940

### Belluscura, Ltd.

Oxford Center for Innovation New Road, Oxford OX1 1BY United Kingdom tel: +44 (0) 1865 261506 toll free: +1 844-883-0940

Copyright © 2017 - Belluscura, Ltd 2017/02 6400-090-001 Rev B