



The Belluscura Bariatric Evacuation Slyde provides a safe, efficient way to evacuate non-ambulatory patients from any multistory building using the stairways. It folds up and stores in stacks up to five units in a storage sleeve. The sleeve can be wall-mounted for quick access in key locations.

## ക്ര

## belluscura

- Folds flat and stacks for easy storage
- Restraints with buckles to secure the occupant for transport
- Corrugated polypropylene, double-wall construction is lightweight and durable
- Strong web handles at corners and sides provide multi-point control for two to seven carriers
- The length of the web handle allows an ergonomically correct upright posture for operators during stair descent

**800LB** Capacity



Copyright © 2017 - Belluscura plc 2017/02 6400-090-002 Rev C



#### **Features and Benefits**

- Folds flat and stacks 5 deep for easy storage
- Easy set-up and instruction pictograms
  on each unit
- Four strong web handles for multipoint control by 2 to 7 operators
- Size permits continued stairwell traffic flow during occupant evacuation

- Functions as post-evacuation surge capacity device
- Lightweight for easier transports
- Buckle restraints are quick and easy to use<sup>1</sup>
- Compatible with Belluscura Stairwell Belay System for heavy occupant evacuations.
- Wall-mounted storage options are space-saving and convenient





Model Number Material Height Width2 Length Weight Folded Height Max. Weight Capacity3,4 6410-000-000 Corrugated polypropylene 8 in. (20.3 cm) 24 in. (61 cm) 85 in. (215.9 cm) 9 lb. (4.1 kg) 2 in. (5.1 cm) 800 lb. (363 kg)

### ഹ്ര

# belluscura

- 1. Passenger restraint set includes waist and chest restraints.
- 2. Dimensions are measured from the outermost edges of the main frame. Specifications are rounded to the nearest whole number. Conversions are calculated before rounding.
- 3. Operators must consider the weight of the passenger, equipment and accessories when determining the total weight on the product.
- 4. Maximum weight capacity is total weight distributed in accordance with basic human anatomy. Belluscura reserves the right to change specifications without notice.

### **ORDER INFORMATION**

email: customerservice@belluscura.com tel: +1 214-396-8313 toll free: +1 844-883-0940

**Belluscura LLC** 

5504 Democracy Drive, Suite 200 Plano, TX 75024 United S tates tel: +1 214-396-8313

#### **Belluscura plc**

82 St John Street, Farringdon London EC1M 4JN United Kingdom tel: +44 (0) 1865 261506 toll free: +1 844-883-0940

Copyright © 2017 - Belluscura plc 2017/02 6400-090-002 Rev C